

PRODUCT SPECIFICATION

REF 51011010 SOFT FOAM BANDAGE 4,5 M



1 GENERAL PRODUCT DESCRIPTION

Description

Cederroth Soft Foam Bandage is an elastic, self-adhesive, latex free bandage made of polyurethane foam. The bandage is protected by paper on both sides and packed in a cardboard box with dispenser function.

The bandage is intended to be used on minor wounds. The skin friendly, breathable bandage does not adhere to the wound nor skin and it does not lose its adhesive properties in water. It can also be used together with ointments or cream.

The bandage is clearly visible thanks to its blue colour which makes it suitable for use e.g. in the food industry.



The product is intended for single use.

CE labelled Medical device, class I.

Storage information

Dry storage in room temperature. For maximum protection, keep the products in their original packages until use.

2 DETAILED PRODUCT INFORMATION

Content per consumer unit

Product/Material

Soft Foam bandage packed in cardboard box

Quantity

4,5m

Dimensions (w x l mm)

Bandage: 60 x 4500

Material

Bandage: Blue polyurethane foam without latex.

Box: Printed cardboard.

Shelf life: 5 years from production date

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310610017293	7310800017454	7312020018112
Content:	1 roll	10 CU	10 SKU/layer x 8 layers
Size (w x l x h mm):	118 x 65 x 108	335 x 243 x 123	800 x 1200 x 1134
Gross weight:	52 g	676 g	79 kg
Packing material:	Carton 17 g	Corrugated board 156 g	-

