

PRODUCT SPECIFICATION

REF 51011003 CEDERROTH FIRST AID & BURN STATION



1 GENERAL PRODUCT DESCRIPTION

Description

A First Aid station containing a variety of first aid products, intended to be used in first aid situations, especially suitable for work places where there is a risk of burns. All products are easily available and protected against dust and moisture by transparent lids. Space for additional products is provided behind a lid with first aid instructions. Adhesive bandages (plasters) and wound cleansers are placed in dispensers for easy access and theft-proof storage. A special key is provided for refill replacement.

Medical devices in the station are CE-labelled.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product</u>	<u>Qty</u>	<u>REF</u> <u>(refill)</u>	<u>Product</u>	<u>Qty</u>	<u>REF</u> <u>(refill)</u>
First Aid station (empty)	1	-	Burn Gel Dressing 2-p	2	901900
Cederroth 4-in-1	2	1910	Burn Gel Spray	1	51011005
Cederroth 4-in-1 Mini	3	1911	Salvequick Wound Cleanser 20-p	1	323700
Salvequick Blue Detectable Plasters, regular, in dispenser	1	51030127	Cederroth Soft Foam Bandage, blue, in dispenser	1	51011011
Salvequick Blue Detectable Plasters, regular/fingertip, in dispenser	1	51030126	First Aid Instruction (label)	1	-
Cederroth Breathing mask	1	1921	Key to the dispensers	1	51616620

Material

First aid station:	Green parts:	Polystyrene, ABS-plastic
	Transparent lid:	Polycarbonate
	Sealing profiles:	EPDM rubber (synthetic rubber)
First aid products:	Refer to Product specifications for each product	

Shelf life: Refer to Product specifications for each product in the station.

3 LOGISTIC INFORMATION

	STATION	SKU	Pallet
GTIN:	-	7310800015412	7312020016026
Content:	Multiple products, see above	1 Station	11 SKU/layer x 3 layers
Size (w x l x h mm):	286 x 117 x 558	294 x 572 x 123	800 x 1200 x 1032
Gross weight:	2925 g (full station)	3256 g	132 kg
Packing material:		Carton 313 g Plastic 18 g	-

