



PRODUCT SPECIFICATION



Product: SALVEQUICK WOUND CLEANSER
Reference code: 3239



1. General description

Sterile sachets made of non-woven tissues, which are saturated in a wound cleansing liquid. The tissues are packed in a well-sealed foil, which prevents leakage and guarantees sterility.

2. Intended use

Salvequick Wound Cleanser is intended for mechanical cleansing of minor wounds when clean tap water is not available.

3. Technical description

Product dimensions

Tissue: 120 x 110 mm
Sachet: 60 x 58 mm

Material

Tissue: Nonwoven made of viscose
Liquid: Physiological saline solution
Foil: Paper/Aluminum/Surlyn® laminate

4. Sterilization

Sterilization method: Irradiation

5. Composition

Transport package

Number of sachets in transport package: 1000
Dimensions (l x w x h) cm: 39 x 29 x 170
Weight (kg): 4,2

6. Labeling

Product

Cederroth contact information
Product name, Reference code, Lot-number,
Expiry date
Symbols: CE, DO NOT REUSE, LOT, USE BY

Transport package

Cederroth contact information,
Product name, Reference code, Lot-
number, Expiry date, EAN-code,
Quantity.

7. Storage conditions

The product should be stored in the original package for maximum protection until time of use.

Shelf life: 5 years from production date

Authorization: