

PRODUCT SPECIFICATION

REF 323700 SALVEQUICK WOUND CLEANSER 20-PACK



1 GENERAL PRODUCT DESCRIPTION

Description

Sterile sachets made of non-woven tissues, which are saturated in a wound cleansing liquid. The tissues are packed in a well-sealed foil, which prevents leakage and guarantees sterility. Salvequick Wound Cleanser Sachet is intended for mechanical cleansing of minor wounds when clean tap water is not available.

The product is a medical device, class IIa.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

Product

Salvequick Wound Cleanser (Ref# 3239)

Tissue
Sachet

Quantity

20

Dimensions (mm)

100 x 140
60 x 58

Material

Tissue: Nonwoven made of viscose
Liquid: Physiological saline solution. Fill volume: 2,7 ml
Sachet: Paper/LDPE/Aluminium/Surlyn®

Sterilization: Gamma irradiation

Shelf life: 5 years from production date

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310610032272	7310803227010	-
Content:	20 sachets	20 CU	19 SKU / layer / 6 layers
Size (w x l x h mm):	762 x 62 x 61	135 x 130 x 160	Pallet 800 x 1200 x 960
Gross weight:	85 g	1841 g	236 kg
Packing material:	9 g	Carton 130 g	-

