

# PRODUCT SPECIFICATION

REF 3227 SALVEQUICK WOUND CLEANSER 40-PACK



## 1 GENERAL PRODUCT DESCRIPTION

### Description

Sterile sachets made of non-woven tissues, which are saturated in a wound cleansing liquid. The tissues are packed in a well-sealed foil, which prevents leakage and guarantees sterility. Salvequick Wound Cleanser Sachet is intended for mechanical cleansing of minor wounds when clean tap water is not available.

The product is a medical device, class IIa.

### Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



## 2 DETAILED PRODUCT INFORMATION

### Content per consumer unit

#### Product

Salvequick Wound Cleanser (Ref# 3239)

Tissue

Sachet

#### Quantity

40

#### Dimensions (mm)

100 x 140

60 x 58

### Material

Tissue: Nonwoven made of viscose  
Liquid: Physiological saline solution. Fill volume: 2,7 ml  
Sachet: Paper/LDPE/Aluminium/Surlyn®

**Sterilization:** Gamma irradiation

**Shelf life:** 5 years from production date

## 3 LOGISTIC INFORMATION

	CU	SKU	Pallet
<b>GTIN:</b>	7310610032272	7310803227010	-
<b>Content:</b>	40 sachets	10 CU	19 SKU / layer / 6 layers
<b>Size (w x l x h mm):</b>	152 x 62 x 61	135 x 130 x 160	Pallet 800 x 1200 x 960
<b>Gross weight:</b>	169 g	1813 g	230 kg
<b>Packing material:</b>	Carton 15 g	Carton 130 g	-

