



PRODUCT SPECIFICATION

REF 51011014 CEDERROTH BURN GEL DRESSING
/ FACE MASK 30x40 CM

1 GENERAL PRODUCT DESCRIPTION

Description

A foam dressing / face mask impregnated with burn gel for the immediate first aid treatment of minor burns, scalds and sunburns. The gel rapidly cools and soothes the injury, prevents further tissue damage and offers a rapid and ongoing pain relief. The dressing has slits for eyes, nose and mouth that can be opened if used as a face mask.

The dressings are individually packed in sealed, terminally sterilized pouches, that easily can be torn open. The Cederroth Burn Gel Dressing / Face Mask 30x40cm is included in Cederroth First Aid Burn Kit.

CE-labelled Medical Device.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

Product/Material

Burn gel dressing in sealed pouch, sterile. Single use.

Dimensions (mm)

Dressing 300 x 400

Pouch 230 x 320

Quantity

1

1

Material

Burn Gel: Water (90%) based gel with propylene glycol, glycerol, emulsifier and surfactant.

Dressing: Polyurethane foam, 3 mm thick

Pouch: Laminate of PE/ALU/PET

Box: N/A

Shelf life: 5 years

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310610021818	7310800021987	7312020022621
Content:	1 dressing	3 CU	12 SKU/layer x 13 layers
Size (w x l x h mm):	230 x 320 x 16	188 x 338 x 80	Pallet 800 x 1200 x 1190
Gross weight:	382 g	1276 g	224 kg
Packing material:	Pouch in PE/ALU/PET laminate	Corrugated board 130 g	-

