

PRODUCT SPECIFICATION

REF 51011008 FIRST AID & BURN STATION SET DIN 13157



1 GENERAL PRODUCT DESCRIPTION

Description

A set with a First Aid & Burn station and a First Aid kit containing a variety of first aid products, intended to be used in emergency situations. The first aid products in the kit are according to DIN 13157.

In the station all products are easily available and protected against dust and moisture by transparent lids. Space for additional products is provided behind a lid with first aid instructions.

The DIN 13157 products are easily available inside the well-organized kit. The bright green and yellow coloured kit is clearly visible at a distance. The strong and durable case has a rubber handle which will not slip out of the hand. Easy to bring to an accident at a further distance.

Wall hanger to attach the kit on the wall is included.

Medical devices in the station and kit are CE-labelled.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit:

<u>Product</u>	<u>Qty</u>	<u>REF (refill)</u>
Cederroth First Aid & Burn Station:	1	51011003
DIN Kit:	1	51611220
Contents:		Refer to Product specifications for each product

Material:

First aid station:	Green parts:	Polystyrene, ABS-plastic
	Transparent lid:	Polycarbonate
	Sealing profiles:	EPDM rubber (synthetic rubber)
DIN Kit:	Case:	PU coated EVA with 75D fabric inside.
	Inner pockets:	Polyethylene
	Handle:	PVC (phthalate-free)
	Zipper:	300 D Nylon and plastic, double puller made of zinc alloy. PES threads.
	Wall hanger:	Green metallic
Contents:		Refer to Product specifications for each product
Shelf life:		Refer to Product specifications for each product in the station.

3 LOGISTIC INFORMATION

	STATION & KIT	SKU	Pallet
GTIN:	Not used	7310800017430	7312020018099
Content:	1 station and 1 kit	1 set of station and kit	4 SKU/layer / 6 layers
Size (w x l x h mm):	Station: 286 x 117 x 558 Kit: 310 x 80 x 247	315 x 572 x 227 (123 + 104)	800 x 1200 x 1200
Gross weight:	Station: 2943 g (full station) Kit: 1562 g (full kit incl. wall hanger)	4977 g	144 kg
Packing material:	-	Carton 462 g Plastic 10 g	-



Version: 01
Issued: 2019-01-11