

Product: STRAP FOR CEDERROTH WALL BRACKET FOR EYE WASH

Reference code: 7205



1. General description

A black strap with Velcro® fastening.

Accessory part to the Cederroth Wall Bracket for Eye Wash (REF 7200).

The strap has a red grip for fast and easy opening of the Velcro® fastening

The strap should be fixed between the Wall Bracket and the surface behind before the Wall Bracket is screwed onto the surface.

2. Intended use

Intended to hold the Cederroth Eye Wash (REF 7251) fixed in place in the Cederroth Wall Bracket (REF 7200) when used in vehicles like trucks, vans, lorries, timber machines and other mobile work places.

3. Technical description

Product dimensions : Strap: (w x l) mm: 20 x 355

Material: Strap: Strong , woven polypropylen

4. Composition

Consumer package

Transport package

<u>Ref#</u>	<u>Product/Material</u>	<u>Pcs</u>	Number of consumer packages
51618205	Strap	7	in transport package: 10
51618662	Plastic bag	1	Dimensions (l x w x h) cm: 21 x 15,5 x 11
			Weight (kg): 0,18

5. Labeling

On product/ consumer package

Cederroth address, telephone- and faxnumber, web-address
Product name, Reference code
EAN-code

On transport package

Cederroth address, telephone- and faxnumber, web-address, Product name, Reference code
Lot-number, EAN-code, Quantity,
Symbols: LOT

6. Storage conditions

The product should be stored in the original package for maximum protection until time of use.

7. Accessories

The strap is intended to be used together with the wall bracket

Authorization: