



PRODUCT SPECIFICATION



Product: SALVEQUICK
WOUND CLEANSERS AND PLASTERS
Reference code: 1880



1. General description

Wound Cleansers and adhesive bandages (plasters) packed in a plastic zip bag.

2. Intended use

To be used in first aid situations. The Wound Cleansers are intended for mechanical cleansing of minor wounds. The adhesive bandages are intended to be used as mechanical barrier and for absorption of exudates from minor wounds.

3. Technical description

Dimensions (mm): Product: 120 x 180 x 20 Adhesive bandages: 25 x 72

Material

Wound cleansers	Sterile tissues soaked in 0.9% NaCl. <i>(for more info, refer to separate product specification REF 3239)</i>
Adhesive bandages	Microperforated polyethylene foil coated with polyacrylate adhesive. Non-woven wound pad (viscose and PP/PE) with a non-stick wound contact layer (PE). Wrapped in cold seal paper.
Zip bag	Polyethylene

4. Composition

Consumer package

Ref#	Product/Material
3239	Wound Cleansers
984213	Adhesive bandages senior
51618261	Plastic zip bag

Transport package

Pcs	Consumer packages in transport package:	
10		20
20	Dimensions (l x w x h) cm:	26 x 21 x 18
1	Weight (kg):	1,2

5. Labeling

Consumer package

Cederroth contact information, Product name, Reference code, CE, Lot-number, Expiry date, EAN-code, Quantity

Transport package

Cederroth contact information, Product name, Reference code, CE, Lot-number, Expiry date, EAN-code, Quantity

6. Storage conditions

The product should be stored in the original package for maximum protection until time of use.
Shelf life: 5 years from production date for each product

Authorization: