

PRODUCT SPECIFICATION

REF 390100 FIRST AID KIT SMALL



1 GENERAL PRODUCT DESCRIPTION

Description

A small First Aid kit for personal use containing first aid products to be used in emergency situations (bleeding). The convenient, soft cover can be attached to a belt or carried in a pocket. All products are easily available and protected against dust, dirt and moisture. The bright green and yellow coloured kit is clearly visible at a distance.

The Medical device products inside are CE-labelled.

Storage information

Dry storage in room temperature. For maximum protection, keep the products in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product/Material</u>	<u>Quantity</u>
Cover	1
Salvequick Plasters (72 x 25 mm)	6
Salvequick Wound cleansers, sterile	2
Salvequick Maxi Cover plaster (76 x 54 mm)	2
Cederroth 4-in-1 Bloodstopper	1

Material

Cover: Nylon 600D Fabric with PU coating on inside
First Aid products: Refer to Product specifications for each product

Shelf life: 5 years from production date for each component

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310613901001	7310803901002	7312023901008
Content:	1 kit	12 CU	21 SKU/layer x 6 layers
Size (w x l x h mm):	145 x 51 x 74	260 x 166 x 168	Pallet 800 x 1200 x 1158
Gross weight:	84 g	1119 g	166 kg
Packing material:	Carton 8 g	Corrugated board 111 g	-

