

PRODUCT SPECIFICATION

REF 51011015 CEDERROTH BURN GEL DRESSING
20x20 CM



1 GENERAL PRODUCT DESCRIPTION

Description

A foam dressing impregnated with burn gel for the immediate first aid treatment of minor burns, scalds and sunburns. The gel rapidly cools and soothes the injury, prevents further tissue damage and offers a rapid and ongoing pain relief.

The dressings are individually packed in sealed, terminally sterilized pouches, that easily can be torn open. The Cederroth Burn Gel Dressing 20x20cm is included in Cederroth First Aid Burn Kit.

CE-labelled Medical Device.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product/Material</u>	<u>Dimensions (mm)</u>	<u>Quantity</u>
Burn gel dressing in sealed pouch, sterile. Single use.	Dressing 200 x200	1
	Pouch 150 x 180	1

Material

Burn Gel:	Water (90%) based gel with propylene glycol, glycerol, emulsifier and surfactant.
Dressing:	Polyurethane foam, 3 mm thick
Pouch:	Laminate of PE/ALU/PET
Box:	N/A

Shelf life: 5 years

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310610021825	7310800021994	7312020022638
Content:	1 dressing	6 CU	35 SKU/layer x 6 layers
Size (w x l x h mm):	150 x 180 x 17	155 x 210 x 105	Pallet 800 x 1200 x 1080
Gross weight:	160 g	1035 g	242 kg
Packing material:	Pouch in PE/ALU/PET laminate	Corrugated board 75 g	-

